


Peak shape issues at the apex (top) of a peak are typically related to detector saturation. This symptom can be observed as flattening and/or noise at the peak apex.



## Preparative Scales

Peak disturbances are commonly observed in preparative separations. Such a symptom is often considered permissible when the chromatogram is used only as a time marker for fraction collecting, rather than quantification purposes.

## Analytical Problems

For quantitative purposes, this will be a problem as the peak area/height will fall above the linear calibration range of the detector.

Injecting a lower concentration within the calibration range will be necessary. If subsequent

impurities present in low amounts cannot be quantified at their required levels with a lower injection concentration, a “high-load:low-load” injection scheme could be used.

With a “high-load:low-load” injection scheme, the same sample is injected at the high concentration for the impurity quantification and at the low concentration for assay level quantification.

## Beer-Lambert Law

The Beer-Lambert Law relates the absorption of light to the properties of material through which the light is traversing.

For UV/Vis detection in chromatography, it is often simplified to  $A = \epsilon cl$ , where  $A$  is the absorbance (inversely proportional to the transmittance that is measured),  $\epsilon$  is the extinction co-efficient (for a particular molecule eluted in a given set of conditions),  $c$  is the concentration of the molecule(s), and  $l$  is the path length through which the light is traveling.

The extinction co-efficient and path lengths are therefore considered constant, leaving the absorbance directly proportional to the concentration.

## Troubleshooting

If this is observed under method conditions previously established without the problem, likely contributing causes could be an incorrect dilution factor in the sample preparation, a different detector with a lower range being used, or an otherwise higher than expected sample concentration in the detection path.

As seen in the Beer-Lambert equation, changing the detection path length with a different

sized flow cell could also change the absolute responses.

---

**Related resources:**

- [UHPLC/HPLC Columns](#)
- [HPLC Column Match Web Tool](#)
- [HPLC Animation Video](#)

Share with friends and coworkers:

- [Click to email this to a friend \(Opens in new window\)](#)
- [Click to share on Twitter \(Opens in new window\)](#)
- [Click to share on Facebook \(Opens in new window\)](#)
- [Click to share on Pinterest \(Opens in new window\)](#)
- [Click to share on LinkedIn \(Opens in new window\)](#)
- [Click to share on Tumblr \(Opens in new window\)](#)
- [Click to share on Reddit \(Opens in new window\)](#)