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Conditions

Column: RP C8 5µm, 4.6 x 150mm

Elution: 0 – 60% ACN with 0.1%TFA over 60 min at 1 ml/min.

Temperature: As indicated.

Sample: Tryptic digest of human growth hormone.

Reference: Hancock, Chloupek, Kirkland and Snyder, J. Chromatography 686, 31-43 (1994)

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APPLICATIONS

Effect of Temperature on Single Stranded Oligonucleotide Analysis

Ivan Lebedev, Dani Xing, Brian Rivera
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501, USA

Overview

Therapeutic oligonucleotides represent a recent breakthrough in the pharmaceutical industry. However, characterization of oligos, specifically by ion-pair reversed phase liquid chromatography (IP-RPLC), can be quite challenging. Oligos are manufactured by solid phase synthesis, where nucleotides are added in a step-wise manner. As such, impurities such as n-1 and n+1 must be characterized, and this may require extensive method development to optimize. Further, it may be necessary for characterization and quantitation of other closely related impurities, such as oxidation of phosphorothioates.

Another challenge is the variation in chromatography seen with oligo analysis. Retention times, peak shapes, and peak area recovery are often variable from injection to injection in any given sequence. One assumed cause of these variations is intramolecular interactions that compromise chromatographic performance. As such, high temperatures, exceeding 60 °C, are often implemented to improve separation. Although this is a requirement when working with double-stranded oligos such as siRNAs, there is still a question of whether this is also necessary for single stranded oligos. Here we present the effect of temperature on two oligos, a 5' conjugated oligo and a phosphorothioate, and how it might be implemented for method development of single stranded oligonucleotides.

Figure 1 illustrates the effect of increasing the method running temperature in 5 °C increments. Retention time for the 5'-Amino C12 oligonucleotide decreases as temperature increases, which is common for any chromatographic method. Often, efficiency and peak shapes are improved at higher temperatures, though that is not observed in this example. Further, in comparing the 45 °C and 65 °C impurity profiles as shown in the insets show little to no differences. One might conclude then that selectivity is not being effected.

To better confirm selectivity changes, mass spectrometry is necessary to identify and characterize impurity peaks. A 22 mer DNA Phosphorothioate, with an unknown sequence was run by high resolution MS. The measured mass for the thioate was 6772.6 Da. The oligo was then run at different temperatures to observe any changes in selectivity for the impurity profiling.

In comparing methods run at 60 and 70 °C (**Figure 2**), as expected, there was a marked decrease in retention for the main peak with the higher temperature running conditions. However, both methods gave similar impurity profiles, with the three earlier eluting impurities observed. The deconvoluted mass spectrum confirmed the same elution order; Impurity Peak 1 was 6443.4 Da, Peak 2 was 6468.4 Da, and Peak 3 was 6170.8 Da. Deconvoluted spectra for Peak 1 with both running conditions are shown in **Figure 3**. Interestingly, Peak 1 and 2 are likely both n-1 impurities, with Peak 1 being target sequence minus guanosine, and Peak 2 being target sequence minus thymidine.

In summary, temperature is often utilized for oligonucleotide analysis. However, for single stranded oligos, there may not be a benefit to running at temperatures exceeding 60 °, as increases in temperature may not improve chromatography nor effect selectivity as one might observe with other macromolecules.

LC Conditions

- Columns:** bioZen™ 2.6 µm Oligo
Dimension: 100 x 2.1 mm (Figure 1) [00D-4790-AN](#)
150 x 2.1 mm (Figures 2,3) [00F-4790-AN](#)
Mobile Phase: Figure 1:
A: 12.5 mM HFIP, 4 mM TEA in Water
B: 12.5 mM HFIP, 4 mM TEA in Methanol
Figures 2,3:
A: 100 mM HFIP, 4 mM TEA in Water
B: 100 mM HFIP, 4 mM TEA in Methanol
Gradient: 5-30 % B in 14 minutes (Figure 1)
25-95 % B in 15 minutes (Figures 2,3)
Flow Rate: 0.3 mL/min
Injection: 1 µL
Temperature: As noted in Figures
Detection: UV @ 260 nm (Figure 1,2)
TOF-MS (Figure 3)
Sample: 5'-Amino C12 Oligo (Figure 1)
DNA 22 mer Phosphorothioate (Figure 2,3)

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
Luna Omega Stability

Mobile Phase	Temperature Limit
0.1% Formic acid	90° C
0.1% TFA	80° C
Phosphate at pH 6	50° C
Phosphate at pH 7.5	50° C/30° C*
Ammonium Bicarbonate pH 8.5	50° C/30° C*

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
Share with friends and coworkers:

- [Click to email this to a friend \(Opens in new window\)](#)

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- [Click to share on Twitter \(Opens in new window\)](#)
- [Click to share on Facebook \(Opens in new window\)](#)
- [Click to share on Pinterest \(Opens in new window\)](#)
- [Click to share on LinkedIn \(Opens in new window\)](#)
- [Click to share on Tumblr \(Opens in new window\)](#)
- [Click to share on Reddit \(Opens in new window\)](#)