[Descargar PDF](#)

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
**Figures 5 – 7** compare the trap performance of MicroTrap C18, MicroTrap Polar, and MicroTrap PS in combination with the Kinetex® 2.6 µm XB-C18 and under the same chromatographic conditions.

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## APPLICATIONS

### Comparison of Three Unique and Complementary Micro LC Columns and Three Trap Selectivities Under Reversed Phase LC-MS/MS Conditions

Roxana Eggleston-Rangel, Ryan Splitstone, and Jason Anspach  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

#### Overview

Micro LC uses columns with internal diameters (I.D.) that are typically within the range of 0.3 - 0.5-mm and packed with the same chromatographic materials used in traditional analytical column dimensions. As is the case in analytical scale LC separations, the chromatographic media impacts the chromatographic performance and selectivity of the separation and can be optimized. A common practice in micro LC is to use a trap-and-elute injection mode. The use of a trap-and-elute injection mode allows for both the significant decrease in injection time and the protection of the column from sample contamination. Micro LC traps can influence the recovery and analytes from the sample based upon the stationary phase chemistry chosen.

In this application, we investigate the influence of different micro LC column and trap stationary phase selectivity combinations on chromatographic performance, sample recovery, and selectivity. We compared micro LC columns in a 50 x 0.3 mm column dimension packed with three different micro LC stationary phase chemistries: Kinetex<sup>®</sup> XB-C18, Luna<sup>®</sup> Omega PS C18, and Luna Omega Polar C18.

The comparison was generated using a sample of 20 stable-isotope-labeled (SIL) peptides under general reversed-phase mobile phase conditions and MS/MS detection using a SCIEX<sup>®</sup> 5500 QTRAP<sup>®</sup>. This application highlights the differences in chromatographic performance and recovery when different combinations of micro LC column and trap selectivity are combined in trap-and-elute mode. Adjusting the column and trap selectivity combination is a useful method development tool for the optimization of your micro LC separation.

**Kinetex 2.6  $\mu$ m XB-C18:** A C18 that is modified with protective iso-butyl side chains for improved analysis of polar compounds.

**Luna Omega PS C18:** The surface contains a positive charged ligand which improves peak shape for basic compounds, while the C18 ligand promotes general reversed phase hydrophobic retention.

**Luna Omega Polar C18:** The C18 ligand provides general hydrophobic interactions while a polar modified particle surface provides enhanced polar compound retention

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Share with friends and coworkers:

- [Click to share on LinkedIn \(Opens in new window\)](#)
- [Click to share on Facebook \(Opens in new window\)](#)
- [Click to share on Twitter \(Opens in new window\)](#)
- [Click to share on WhatsApp \(Opens in new window\)](#)
- [Click to email a link to a friend \(Opens in new window\)](#)